



Risk assessment for medical devices

Motivation & Context

Medical devices domain characteristics


- Significant documentation effort to gain certification
- Today, document driven development approach
- Risk management (RM) activities are a cost driver
- Lack of standardized methods and practices to qualify medical equipment
- RM information-wise decoupled from engineering

Objectives

- Provide standard based method implementation for RM
- Enable integration of RM into engineering activities
- Comply with RAAML language and philosophy
- Enable MBSE

Related standards


- ISO 14971:2019 – medical devices – Application of risk management to medical devices
- IEC 62304:2006 Medical devices software – Software lifecycle processes
- IEC 62366:2015 Medical devices - Part 1: Application of usability engineering to medical devices
- IEC 60601:2006 - Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
- RAAML Version 1.0
- SysML Version 1.6
- UML Version 2.5.1



Example


Medicussy in action

Define the system boundary




- The blue colored blocks are specializations of Medicusssy library items



Usability analysis - normal use




Usability analysis - normal use



Usability analysis - use errors



Usability analysis - abnormal use





Abnormal use is defined as
“conscious, deliberate act [..] that
is counter to or violates NORMAL
USE

Source: ISO 62366:2015

Risk modeling

ibd [Situation] Reservoir empty - indication overlooked #1 [Reservoir empty - indication overlooked #1]




Risk analysis matrix

#	△ Name	Event(s)	Hazardous Situation (HS)	Initial HS Probability	Hazard	Causality Potential	Initial Harm Probability	
1	Deliberate insulin overdose	Repeating manual injection command	Insulin overdose		Hypoglycaemic patient			Death
2	Electrostatic device damage - decreased consciousness	Electrostatically charged patient touches device	ESD causes pump and alarms are failir	5.0			5.0	Decrease
3	Electrostatic device damage - organ damage	Electrostatically charged patient touches device ESD causes damage of alarm and pump	ESD causes pump and alarms are failir	5.0			5.0	Minor c
4	Reservoir empty - empty ampoule inserted	Empty reservoir removed Insert empty ampoule Indication of low reservoir	Insulin injection failes		Hyperglycaemic patient			Decrease
5	Reservoir empty - indication overlooked #1	Indication of low reservoir not perceived by patient Indication of low reservoir	Overlook insulin reservoir emptyness	2.0	Hyperglycaemic patient	Patient ignores body signals for <=30 Min	3.0	Decrease
6	Reservoir empty - indication overlooked #2	Indication of low reservoir not perceived by patient Indication of low reservoir	Overlook insulin reservoir emptyness	2.0		Patient ignores body signals for >30 Minut	2.0	Minor c


	Hazardous Situation (HS)	Initial HS Probability	Hazard	Causality Potential	Initial Harm Probability	Harm	Initial Harm Severity	Requirement relations	Mitigation	Risk verification
	Insulin overdose		Hypoglycaemic patient			Death		F 86 Manual injection	F 78 Avoid overdose	
es device	ESD causes pump and alarms are failir	5.0			5.0	Decreased consciousness	3.0		R 73 Electrostatic discharge F 74 Alarm health monitor F 75 Pump health monitor	
es device ip	ESD causes pump and alarms are failir	5.0			5.0	Minor organ damage	5.0		R 73 Electrostatic discharge F 74 Alarm health monitor F 75 Pump health monitor	
	Insulin injection failes		Hyperglycaemic patient			Decreased consciousness	3.0	D 84 Exchangeable ampoule	F 83 Detect empty ampoule	
ed by patient	Overlook insulin reservoir emptyness	2.0	Hyperglycaemic patient	Patient ignores body signals for <=30 Min	3.0	Decreased consciousness	3.0	F 79 Indicate filling level F 80 Alert low reservoir	Ph 82 Vibration actuator F 81 Increase alarm inter	
ed by patient	Overlook insulin reservoir emptyness	2.0		Patient ignores body signals for >30 Minut	2.0	Minor organ damage	5.0			

ISO 14971 Profile & Library

Medicussy structure





ISO 14971 Profile



MedicalRiskItem is not a normative term within ISO 14971.


The ISO 14971 profile introduces the stereotype MedicalRiskItem in order to manage risk in the model analogue to RAAML FMEA implementation.

ISO 14971 Library



Source: ISO 14971:2019 Annex C

Open issues




General questions

- Recommendations how to share with the community?
- Do you see this project as relevant future RAAML extension?
- Smart modeling with library driven approach, especially:
 - Creating risk model compositions
 - Automatic generalization
 - Redefine assistance

Implementation


Potential addons to support risk management method implementation

Medicussy Library




Selected library items that are introduced to support medical system analysis in accordance to applicable ISO 14971

Medicussy Profile




Selected stereotypes that are introduced to support medical system analysis

IEC 62304 & IEC 62366 Library



Selected library items that are introduced to support medical system analysis in accordance to further applicable standards



Risk management process – initial analysis

